

# B E D R O C K

## B I O S C I E N C E

### **PRO SERIES** **FOOT LIGHT**

The Bedrock Bioscience Light Pad is a flexible infrared light device designed to increase circulation, decrease pain, and promote healing.

## QUICK START GUIDE

Plug in and press both power buttons to turn on lights and magnet.  
DO NOT STAND on lights.

## FOR STORAGE

**DO NOT FOLD LIGHTS, ONLY ROLL** as damage will occur to the light unit and void the warranty.  
Store in the carrying case properly when not in use.

## CLEANING CARE

Use an alcohol wipe or spray on cloth and wipe down.  
DO NOT spray directly onto unit or wash in water.