



B E D R C K

B I O S C I E N C E

L I G H T P A D

The Bedrock Bioscience Light Pad is a flexible infrared light device designed to increase circulation, decrease pain and promote healing.

QUICK START GUIDE

Plug in and hold power for 3 seconds to turn on lights.
During treatment, DO NOT STAND on lights.

FOR STORAGE

DO NOT FOLD LIGHTS, ONLY ROLL as damage will occur to the light unit and void the warranty.

CLEANING CARE

Use an alcohol wipe or spray on cloth and wipe down.
DO NOT spray directly onto unit or wash in water.