

# B E D R C K

## B I O S C I E N C E

### BODY REBOOT LIGHT P A D

The Bedrock Bioscience Body Reboot Light Pad is a flexible LED light device designed to increase circulation, decrease pain, and promote recovery.

### QUICK START GUIDE

Plug in and hold power for 3 seconds to turn on lights.  
DO NOT STAND on lights.

### FOR STORAGE

**DO NOT FOLD LIGHTS, ONLY ROLL** as damage will occur to the light unit and void the warranty.

### CLEANING CARE

Use an alcohol wipe or spray on cloth and wipe down.  
DO NOT spray directly onto unit or wash in water.